

# Gel Shift/ EMSA Protocol

**Jianping Ye, MD**  
**Antioxidant and Gene Regulation Lab**  
**Pennington Biomedical Research Center**  
**Louisiana State University**  
**6400 Perkins Road**

## References:

1. **Ye J**, Ghosh P, Cippitelli M, Subleski J, Hardy KJ, Ortaldo JR, Young H. Characterization of a silencer regulatory element in the IFN- $\gamma$  promoter. *Journal of Biological Chemistry* 1994;269(41):25728-25734.
2. **Ye J**, Zhang X, Dong Z. Characterization of the human GM-CSF promoter: AP1 and a Sp1-related protein activate the promoter activity that is suppressed by YY1. *Molecular and Cellular Biology* 1996;16(1):157-167.
3. **Ye J**, Cippitelli M, Dorman L, Ortaldo JR, Young HA. The nuclear factor YY1 suppresses the human IFN- $\gamma$  promoter through two mechanisms: inhibition of AP1 binding and activation of a silencer element. *Molecular and Cellular Biology* 1996;16(9):4744-4753.

## Electrophoresis System

Vertical gel electrophoresis apparatus (BIBCO BRL model V16-2)

### Cast non-denature PAGE Gel:

Stock Solution	50 ml			100 ml		
	3.5%	5%	6%	3.5%	5%	6%
40% acrylamide	4.38ml	6.25ml	7.5ml	8.76ml	12.5ml	15ml
2% Bis-acrylamide	2.91ml	4.15ml	5ml	5.82ml	8.3ml	10ml
10 x TBE buffer	2.5ml	2.5ml	2.5ml	5ml	5ml	5ml
TEMED	50 $\mu$ l	50 $\mu$ l	50 $\mu$ l	100 $\mu$ l	100 $\mu$ l	100 $\mu$ l
10% ammonium persulfate	0.6ml	0.6ml	0.6ml	1.2ml	1.2ml	1.2ml
dH <sub>2</sub> O	39.6ml	36.5ml	34.4ml	79.1ml	72.9ml	68.7ml

### Prepare 0.5 x TBE running buffer 1800ml:

10 x TBE buffer      90ml  
dH<sub>2</sub>O                    1710ml

#### 10 x TBE formula:

Tris base                    108 g (89 mM)  
Boric acid                    55 g (89 mM)  
0.5 M EDTA (pH 8.0)      40 ml  
Use dH<sub>2</sub>O to bring total volume to 1000 ml

### Prerun the gel:

Run the gel at 200 V for 30 min

(Prepared by Jianping Ye, M.D.)

### **Label DNA probe with T4 kinase:**

DNA probe	5 $\mu$ l (100 - 500 ng)
5 x Forwarding buffer	4 $\mu$ l
<sup>32</sup> P- $\gamma$ -ATP	5 $\mu$ l
dH <sub>2</sub> O	5 $\mu$ l
T4 kinase (10 Unit/ $\mu$ l)	1 $\mu$ l

Incubate at 37 °C for 30 min

### **Label DNA probe with Ready-To-Go kit (Pharmacia Biotech: Cat# 27-5335-01):**

1. Use dH<sub>2</sub>O to dissolve reagent: 40  $\mu$ l/vial
2. Add in DNA probe (100 - 500 ng): 5  $\mu$ l
3. Add in <sup>32</sup>P- $\alpha$ -dCTP 5  $\mu$ l
4. Incubate at 37 °C for 15 min

### **Purification of labeled probe with G50 micro column (Pharmacia Biotech 27-5335-01)**

1. Spin 3000 rpm in 5415C eppendorf centrifuge for 1 min to get ride of buffer.
2. Load probe onto the resin
3. Collect purified probe by spinning the loaded column at 3000 rpm in 5415C eppendorf centrifuge for 2 min: use 1.5 ml microcentrifuge tube to collect probe.
4. Check probe activity with scintillation counter: 1  $\mu$ l probe/4ml scintillation solution

### **Prepare the nuclear protein extract:**

1.  $1 - 5 \times 10^7$  cell pellet in 1.5 ml microcentrifuge tube
2. Resuspend cells in 300  $\mu$ l lysis buffer, and keep on ice for 4 min to break the cell membrane
3. Spin cells at 10000 rpm for 1 min at 4 °C to separate nuclei with cytoplasmic component
4. Remove supernatant as cytoplasmic extract
5. Wash the nuclear pellet with 300  $\mu$ l washing buffer
6. Spin cells at 10000 rpm for 1 min at 4 °C to pellet the nuclei
7. Resuspend the nuclear pellet in 100  $\mu$ l (for  $1 \times 10^7$  cells) or 200  $\mu$ l (for  $5 \times 10^7$  cells) extraction buffer
8. Freeze the tube in -70 °C or continue the extraction.
9. Pipette the nuclear resuspension.
10. Spin cells at 14000 rpm for 5 min at 4 °C to collect supernatant as the nuclear extract.
11. Measure protein concentration and adjust it to 1  $\mu$ g/ $\mu$ l with the extraction buffer for use in gel shift assay.

(See the attached formula for each buffer).

M.D.)

## **2 x gel shift reaction buffer**

Stock Solution	15 ml Buffer	Final Concentration
50% Glycerol	3.6 ml	12%
1 M HEPES (pH 7.9)	360 ul	24 mM
1 M Tris-HCl (pH 8.0)	120 ul	8 mM
0.5 M EDTA (pH 8.0)	60 ul	2 mM
100 mM DTT	150 ul	1 mM
dH <sub>2</sub> O	10.7 ml	

## **Prepare reaction mixture for gel shift assay**

1. 2 x reaction buffer 12  $\mu$ l
2. BSA (1  $\mu$ g/ $\mu$ l) 3  $\mu$ l
3. Poly (dI-dC) ( 0.5  $\mu$ g/ $\mu$ l) 2  $\mu$ l
4. Nuclear extract (1  $\mu$ g/ $\mu$ l) 3  $\mu$ l
5. dH<sub>2</sub>O 3  $\mu$ l
6. Keep at Rt or on ice for 10 min without Ab, 20 min with Ab
7. Add in DNA probe (4000 cpm/ $\mu$ l) 1  $\mu$ l
8. Keep at Rt for 20 min
9. Load the gel and run at 200 V for 1 - 1.5 hr. Use DNA loading buffer in lane 1 as indicator of free probe. Free probe usually run at the same mobility as the blue dye of the DNA loading buffer. Stop the gel when the dye runs at 3 cm to the bottom.
10. Dry the gel and expose the dried gel to X-ray film at -70 °C over night.
11. Develop the film.

(Prepared by Jianping Ye, M.D.)

## **Nuclear Protein Preparation Buffer**

### **Lysis Buffer**

<u>Stock Solution</u>	<u>5 ml</u>	<u>10 ml</u>	<u>Final Concentration</u>
1 M KCl	250 ul	500 ul	50 mM
IGEPAL CA-630 (Sigma)	25 ul	50 ul	0.5%
1 M HEPES (pH 7.8)	125 ul	250 ul	25 mM
1 mg/ml Leupeptin (Sigma)	50 ul	100 ul	10 ug/ml
1 mg/ml Aprotinin (Sigma)	100 ul	200 ul	20 ug/ml
250 mM DTT	2.5 ul	5 ul	125 uM
100 mM PMSF	50 ul	100 ul	1 mM
dH <sub>2</sub> O	4.4 ml	8.8 ml	

### **Washing Buffer**

<u>Stock Solution</u>	<u>5 ml</u>	<u>10 ml</u>	<u>Final Concentration</u>
1 M KCl	250 ul	500 ul	50 mM
1 M HEPES (pH 7.8)	125 ul	250 ul	25 mM
1 mg/ml Leupeptin	50 ul	100 ul	10 ug/ml
1 mg/ml Aprotinin	100 ul	200 ul	20 ug/ml
250 mM DTT	2.5 ul	5 ul	125 uM
100 mM PMSF	50 ul	100 ul	1 mM
dH <sub>2</sub> O	4.4 ml	8.8 ml	

### **Extraction Buffer**

<u>Stock Solution</u>	<u>5 ml</u>	<u>10 ml</u>	<u>Final Concentration</u>
1 M KCl	2.5 ml	5 ml	500 mM
1 M HEPES (pH 7.8)	125 ul	250 ul	25 mM
50 % Glycerol	1 ml	2 ml	10%
1 mg/ml Leupeptin	50 ul	100 ul	10 ug/ml
1 mg/ml Aprotinin	100 ul	200 ul	20 ug/ml
250 mM DTT	2.5 ul	5 ul	125 uM
100 mM PMSF	50 ul	100 ul	1 mM
dH <sub>2</sub> O	1.2 ml	2.4 ml	

(Prepared by Jianping Ye, M.D.)

## **Prepare double-stranded DNA probe**

1. Measure OD<sub>260</sub> of single-stranded DNA (5 ul DNA + 995 ul dH<sub>2</sub>O)
2. Calculate the DNA concentration with formula:  
OD<sub>260</sub> x dilution fold (200) x 33 ÷ 1000 = (ug/ul)
3. Adjust DNA concentration to 1 ug/ul with dH<sub>2</sub>O
4. Mix equal amount of the complementary DNA (such as Sense: Antisense = 50 ug: 50 ug)
5. Heat the mixture to 95 °C for 5 min in heating block
6. Leave the heating block at RT to let it cool down slowly
7. Measure OD<sub>260</sub> of the annealed DNA (5 ul DNA + 995 ul dH<sub>2</sub>O)
8. Calculate the DNA concentration with formula:  
OD<sub>260</sub> x dilution fold (200) x 50 ÷ 1000 = ug/ul
9. Dilute some double-stranded DNA to 0.1 ug/ul with dH<sub>2</sub>O
10. Keep double and single-stranded DNA in -20 °C, keep the diluted DNA at 4 °C.

(Prepared by: Jianping Ye, M.D.)